


BloodPorter®

Confidence in your cold chain


User Guide

BloodPorter® Carriers are validated to transport red blood cells and platelets between hospital sites, blood banks and wards, maintaining thermal performance and integrity throughout the journey.

Lid



Product



MediCool® Packs



Inner Chamber



Protective Outer Bag



Prior to use:

- **Red Blood Cells** - pre-condition the MediCool® packs at +5°C for a minimum of 24 hours. Remove 20 mins before use.
- **Platelets** - pre-condition the MediCool® packs at +22°C for a minimum of 24 hours. Remove 20 mins before use.

System assembly and packing

1. Open protective outer bag and inner chamber lid.
2. Place the product and pre-conditioned MediCool® packs into the inner chamber. Ensure the product and cool packs are packed in accordance with the validation (contact Helapet for packing instructions on specific BloodPorter® models).
3. Place the lid of the inner chamber on correctly and seal the protective outer bag immediately.



Please do not open the BloodPorter® Carriers unnecessarily and remember to immediately seal the system. We recommend checking the validation to ensure the carrier is suitable for your intended use.

To re-use the system:

- Conduct a weight, visual and squeeze check of the MediCool® packs to ensure there is no loss of liquid. If the MediCool® is +20g different to the weight shown, we recommend immediate disposal and replacement.
- Chipped or broken inner chambers need to be replaced before use as this can affect thermal performance.



Replacement MediCool® packs and inner chambers are available to purchase separately from Helapet.

Further information is available from:

Helapet

A CliniMed® Group Company

Helapet Limited, Lister House, Blackburn Road, Houghton Regis, Bedfordshire. LU5 5BQ
Telephone 0800 0328 428 Fax +44 (0) 1582 501981
E-mail: sales@helapet.co.uk | Web: helapet.co.uk

