



Sterile Sharpie Marker Pens

Triple-bagged

Our Sterile Sharpie Marker Pens are cleanroom packaged and terminally gamma sterilised for use in cleanroom, isolators and aseptic environments. Triple-bagged polypackaging supports safe and swift aseptic transfer of our pens, reducing the risk of cross-contamination and reducing processing time. Our fine-tip black markers are ideal for the precise labelling of polybags, vials, syringes and labels.

Features and Benefits

- ✓ Permanent writing on most types of surfaces
- ✓ Quick drying, non-toxic black ink
- ✓ Fine tip for precise writing
- ✓ Low-linting, triple-bagged polypackaging supporting aseptic transfer procedures
- ✓ Terminally sterilised for assured use in aseptic environments



Ordering Information

XXTB510BLK
triple-bagged
sterile • 5 x 10 pcs/pck

..... Individually double polybagged



Technical Data

Tip size	Fine
Ink Colour	Black
Properties	Quick drying; resistant to both water and fading; non-toxic
Packaging	10 x individually double-polybagged sharpie pens contained in triple polybagged pack, 5 x triple-polybagged packs per case
Sterilisation	Gamma sterilised
Shelf Life	2 years

Helapet Limited

Lister House | Blackburn Road | Houghton Regis | Bedfordshire | LU5 5BQ

Web: helapet.co.uk | **Phone:** +44 (0)1582 501980 | **Email:** sales@helapet.co.uk

Helapet Limited is a CliniMed Group Company.

